


ATS 40 Transfusion Set

with 40 μ m – Filter



✓ SAFE

✓ EASY

✓ RECOMMENDED

ATS 40 Transfusion Set



- ↪ 40 µm Filter
- ↪ Flexible drip chamber
- ↪ 150 cm length of tubing
- ↪ Easy Roller clamp
- ↪ Luerlock-Connection

Usage: **blood transfusion, specific for autologous transfusion**


CE 0482

List Nr.	Description	units/carton
210-000	ATS-40 Transfusion Set with 40 µm -Filter 150 cm, Roller-clamp, Luer Lock	60

Distribution:

ATS – Autotrans

*Closed Low-Vacuum-System
for postoperative autologous retransfusion*



✓ ***PATIENT SAFE***

✓ ***EASY OF USE***

✓ ***COST EFFECTIVE***

ATS – Autotrans



ATS-Autotrans –
Autotransfusion Set



ATS 40 – Transfusion Set
with 40µ-Filter

- Low-Vacuum – reduced Hemolysis
- Easy of Use
- One closed system
- No device (pump) needed
- Autologous blood reduce usage of homologues red cells
- 1000 ml Retransfusion bag
- 2 One-way-valves
- 200 µm Filter integrated in the system
- Cost effective
- Latex free

CE 0482

List Nr.	Description	units/carton
120-000	ATS-Autotrans Autotransfusion System with 3 Drainage connector and de-connectable Drainage-extension, bellow and 1000ml Retransfusion bag with One-way- valve and 200µ - Filter	10
210-000	ATS-40 Transfusion Set with 40 µm -Filter 150 cm, Roller-clamps, Luer Lock	60

Distribution: