

# AUTOMATIC PRESSURE INFUSOR AUTOPRESS

PRECISE PRESSURE SETTING. PRESET PRESSURE.

BIEGLER **Autopress** is another component in our range of warming equipment. It can be used either directly with our blood and infusion warmers or as a stand-alone solution wherever fluids have to be dispensed under constant pressure.



## AUTOMATIC PRESSURE INFUSOR AUTOPRESS

- Preset pressure is always maintained
- Precise pressure setting due to electronic control
- Accomodated pressure cuffs 2 x 500, 2 x 1000 or 1 x 3000 cc
- Fits on any IV pole
- Small and lightweight
- Pressure limited to 300 mmHg
- Mains or battery operated
- Pressure cuffs can be emptied rapidly and easily
- Significant reduction of set up time
- Battery backup up to 90 minutes

### TECHNICAL DATA:

<b>Device:</b>	automatic pressure infusor
<b>Type designation:</b>	Autopress
<b>Voltage:</b>	9V DC
<b>Power consumption:</b>	max. 12 W
<b>Type:</b>	mains and rechargeable battery operated
<b>Classification:</b>	Ila according to rule 11
<b>Degree of protection against electric shock:</b>	B
<b>Dimension:</b>	90 x 115 x 260 mm
<b>Weight:</b>	1.46 kg
<b>Type of operation:</b>	continuous operation
<b>Mains adaptor Type designation:</b>	12W P2EFMW3
<b>Input:</b>	110-240V AC, 50/60 Hz, 250 mA
<b>Output:</b>	9V DC, 1.33A max.
<b>Accessories:</b>	Pressure cuff 500 cc Pressure cuff 1000 cc Pressure cuff 3000 cc

**AUTOPRESS**

**E. BIEGLER GmbH**  
Allhangstrasse 18a  
A-3001 Mauerbach (Vienna)  
Austria

Tel: +43 1 979 21 05  
Fax: +43 1 979 21 05 16  
office@biegler.com  
www.biegler.com

**BIEGLER**  
MEDIZIN ELEKTRONIK

